

P1 Experiment  
Detection

Ⓢ MW 7:45 pm

Subj. #:

101)

Date: 6/12/97

Subject:

Subj. type: Normal

Amp. Settings (make note of any changes)

	Gain	HPF	LPF	External filter
Ch. 1-28	.10	.01	2000 <del>100</del>	50
Ch. 29	.10	.01	2000 <del>100</del>	
Ch. 30	.10	.01	2000 <del>100</del>	
Ch. 31	.10	.01	2000 <del>100</del>	

Run Calibration

On DIG:

```
cd \dig
digitize sf.hdr <filename>.log <filename>.raw 1
```

On STIM:

```
cd \exprmnts\p1
cal
c 990
s [on DIG]
s [on STIM]
```

After approximately 30 cal pulses,

```
p [on DIG]
q [on STIM]
```

Cap placement

Head circumference: 19 cm

Nasion-to-inion: 37 cm

Channel Assignments: (make note of any irregularities)

1: FPz	9: FC2	17: CP1	25: PO7 +
2: <del>LC</del> (plug-in) <del>RL</del> <del>LE</del>	10: FC6	18: CP2 -	26: PO3 +
3: F3	11: T7 +	19: CP6	27: POz +
+ 4: Fz	12: C3	20: P7	28: PO4
+ 5: F4	13: C4	21: P3 +	29: PO8 +
6: <del>RL</del> <del>LE</del> (plug-in) <del>LC</del>	14: Cz	22: Pz +	30: O1
+ 7: FC5	15: T8 +	23: P4	31: Oz -
+ 8: FC1	16: CP5 +	24: P8	32: O2 +

Comments: