

On STIM

321.CNT

p1 (or p1r or p1rr)

On DIG--Change condition code

	Box #	Comments
Block 1	5 ✓	
Block 2	3 ✓	electrode T8 (15) blocking
Block 3	4 ✓	throughout exp (despite quik-insert)
Block 4	4 ✓	
Block 5	5 ✓	
Block 6	2 ✓	electrode T7 (11) blocking
Block 7	1 ✓	throughout exp.
** break **		
Block 8	1 ✓	
Block 9	5 ✓	
Block 10	2 ✓	
Block 11	5 ✓	
Block 12	3 ✓	
Block 13	4 ✓	
Block 14	3 ✓	
** break **		
Block 15	1 ✓	
Block 16	3 ✓	
Block 17	4 ✓	
Block 18	2 ✓	
Block 19	5 ✓	
Block 20	3 ✓	
Block 21	4 ✓	
** break **		
Block 22	2 ✓	T2 (32) went bad
Block 23	1 ✓	
Block 24	3 ✓	
Block 25	1 ✓	
Block 26	2 ✓	
Block 27	2 ✓	
Block 28	1 ✓	
** break **		
Block 29	5 ✓	
Block 30	4 ✓	
Block 31	5 ✓	
Block 32	2 ✓	
Block 33	3 ✓	
Block 34	1 ✓	
** break **		
Block 35	4 ✓	
Block 36	5 ✓	
Block 37	4 ✓	
Block 38	3 ✓	
Block 39	2 ✓	
Block 40	1 ✓	