

P1 Experiment - *Detection*

~~(Easy Discrimination: O vs. •)~~

Subj. #: 119  
Subject: \_\_\_\_\_  
Subject Category: NC

Date: 8/27/97  
Time: 4:30 pm  
Experimenter: \_\_\_\_\_

**Amp. Settings** (make note of any changes)

|          | Gain | HPF | LPF  | External Filter |
|----------|------|-----|------|-----------------|
| Ch. 1-28 | .10  | .01 | 2000 |                 |
| Ch. 29   | .10  | .01 | 2000 | <u>50</u>       |
| Ch. 30   | .10  | .01 | 2000 |                 |
| Ch. 31   | .10  | .01 | 2000 |                 |

**Run Calibration**

On DIG:

```
cd \dig  
p1 <subj#>
```

On STIM:

```
cd \experiments\p1  
cal  
c 990  
s [on DIG]  
s [on STIM]
```

After approximately 30 cal pulses,

```
p [on DIG]  
q [on STIM]
```

**Cap placement**

Head circumference: \_\_\_\_\_ cm

Nasion-to-inion: \_\_\_\_\_ cm

**Channel Assignments: (make note of any irregularities)**

|                   |         |         |         |
|-------------------|---------|---------|---------|
| 1: FPz            | 9: FC2  | 17: CP1 | 25: PO7 |
| 2: LC (plug-in)   | 10: FC6 | 18: CP2 | 26: PO3 |
| 3: F3             | 11: T7  | 19: CP6 | 27: POz |
| 4: Fz             | 12: C3  | 20: P7  | 28: PO4 |
| 5: F4             | 13: C4  | 21: P3  | 29: PO8 |
| 6: RLOE (plug-in) | 14: Cz  | 22: Pz  | 30: O1  |
| 7: FC5            | 15: T8  | 23: P4  | 31: Oz  |
| 8: FC1            | 16: CP5 | 24: P8  | 32: O2  |

Comments: