

**P1 Experiment  
(Detection)**

Subj. #: 314

Date: 3/30/98

Subject: \_\_\_\_\_

Time: 4:45 pm

Subject Category: FL

Experimenter: \_\_\_\_\_

**Amp. Settings** (make note of any changes)

|          | Gain | HPF | LPF  | External Filter |
|----------|------|-----|------|-----------------|
| Ch. 1-28 | .10  | .01 | 2000 |                 |
| Ch. 29   | .10  | .01 | 2000 | <u>50</u>       |
| Ch. 30   | .10  | .01 | 2000 |                 |
| Ch. 31   | .10  | .01 | 2000 |                 |

**Run Calibration**

On DIG:

```
cd \dig
p1 <subj#>
```

On STIM:

```
cd \experiments\p1
cal
c 990
s [on DIG]
s [on STIM]
```

After approximately 30 cal pulses,

```
p [on DIG]
q [on STIM]
```

**Cap placement**

Head circumference: 51.0 cm

Nasion-to-inion: 31.5 cm

**Channel Assignments: (make note of any irregularities)**

|                        |           |          |           |
|------------------------|-----------|----------|-----------|
| 1: <del>EEG</del> LMAS | + 9: FC2  | +17: CP1 | + 25: PO7 |
| 2: RLOE (plug-in)      | + 10: FC6 | +18: CP2 | + 26: PO3 |
| +3: F3                 | +11: T7   | +19: CP6 | +27: POz  |
| 4: Fz                  | +12: C3   | + 20: P7 | +28: PO4  |
| +5: F4                 | + 13: C4  | + 21: P3 | +29: PO8  |
| 6: LC (plug-in)        | +14: Cz   | + 22: Pz | 30: O1    |
| +7: FC5                | +15: T8   | + 23: P4 | 31: Oz    |
| +8: FC1                | +16: CP5  | 24: P8   | 32: O2    |

Comments:

Lots of alpha!!